



PHS-25 Bench Top PH Tester

(Precision) PHS-25C, 2C, 3C pH meter

(Precision) PHS pH meter is a digital display pH meter, which uses white backlight, which can display pH, or potential (mV) value. The instrument is suitable for laboratory universities, research institutes, environmental monitoring, industrial and mining enterprises and other sectors of water sample to measure pH of the solution and potential (mV) value.

Instrument Model	PHS-25C	PHS-2C	PHS-3C
pH measurement range	0.00 ~14.00pH	0.00 ~14.00pH	0.00 ~14.00pH
pH measurement accuracy	±0.05pH	±0.02pH	±0.01pH
Resolution	0.01pH	0.01pH	0.01pH
mV measuring range	0 ~±1999mV	0 ~±1999mV	0~±1999mV
mV measurement accuracy	±10mV	±10mV	±2mV
Resolution	1mV	1mV	1mV
Temperature compensation range	0~60°C	0~60°C	0~60°C
Standard buffer solution	pH4.00 ; PH6.86 ; pH9.18,		
Power supply type	Power supply adapter, or power cord		
Dimensions	220mm×200mm×85mm		
weight	1KG		